

KBr HYDRAULIC PRESS

'INSTRUMEX' Brand Kbr Hydraulic Press is suitable for making 13 mm dia die OR Ultra micro die with a set of Evocable Kbr Die for IR Spectroscopy solid sampling. The die produces a disc of 13 mm diameter from which a spectrum may be obtained. Using full beam energy of the instrument without use of a beam condenser. A nipple is fitted in the die for evacuated. This die easily can be dismental for cleaning purpose.



ACCESSORIES FOR IR

- * KBr PRESS
- * KBr DIE
- * PLUNGER ANVIL
- * PELLET HOLDER
- * PELLET REMOVER
- * KBr/ NaCl WINDOW/CELL

